




Pro 3D Printers

In-office Dental Manufacturing has Arrived



Fast

Rapidly respond to changing patient needs with flexible in-office manufacturing.

Accurate

Print confidently with repeatable, accurate results backed by university research.

Easy

Spend more time with patients with easy-to-use software and hardware.


Frictionless Workflow

User-friendly, cutting-edge solutions.

SprintRay products are built to work together, creating opportunities for collaboration, innovation, and outstanding clinical results. From design to print to post-processing, SprintRay products help you maximize the efficiency of your clinic.

User-friendly Software

Dashboard 2.0 and RayWare simplify the process of printing and designing appliances. They're easy to use and fully integrated with SprintRay printers.



High Performance 3D Printing

Rapidly print all kinds of dental appliances for treatment with SprintRay Pro. Advanced software and hardware makes 3D printing a breeze.

Powerful Post Processing

Pro Wash/Dry and Pro Cure make a dynamic pair for the final step of your workflow. Get the best out of your parts with SprintRay post processing.



Industry-leading Materials

SprintRay resins offer fantastic value for money. For maximum flexibility, we also certify top industry brands to work with our printers.

More information here.
www.sprinray.com

Get in touch with us.
info.eu@sprinray.com

Pro 95 order No.
SRI-0102004

Pro 55 order No.
SRI-0102014



SprintRay Pro 3D Printer

The Gold Standard for In-office Dental Manufacturing



Intensity, Meet Accuracy

The custom projector in SprintRay Pro delivers high-intensity light evenly across the build surface for perfect, accurate prints every time.

Bespoke 405nm LED projector



Always Cutting-edge

SprintRay Pro features advanced onboard electronics, enabling over-the-air updates that improve the printing experience. Free over-the-air updates bring meaningful improvements.

Free over the air updates



Industrial Throughput

Enough throughput to handle the largest jobs. SprintRay Pro features a large build plate and rapid print speeds to meet the needs of even the most demanding offices.

Up to 3 inches per hour print speed



Open Certified Materials

Our unique resin ecosystem allows you to leverage the industry's greatest innovations, whether they come from us or one of our materials partners.

Over 50 resins to choose from



Pro 55

Precision Yield

Incredible precision makes this unit ideal for all manner of restorative applications.

A 55-micron XY pixel size and Pixel Toning brings industry-leading accuracy to Pro 55, making it perfect for jobs that require razor-thin margins such as temporary and definitive crowns, hybrid dentures, and more. Pro 55 is the ideal printer for clinics or labs that specialize in restorative applications, reducing treatment cost and delivery times while improving your workflow.



Surgical Guides



Hybrid Dentures



Temporary and Definitive Crowns




Cast Applications

Pro 95

Production Workhorse

Accurate production and high throughput for all manner of dental appliances.

SprintRay Pro 95 is the perfect 3D printer for dental offices that want to deliver on the promises of in-office production. Featuring unmatched throughput, high accuracy, and a high flexibility, Pro95 is the right choice for most dental practices. Incredible accuracy, repeatability, and user-friendliness means you can always deliver an outstanding patient experience.



Surgical Guides




Study and
Waxup Models



Models for Aligner
Production



Indirect Bonding
Trays



Digital Dentures
and Try-ins



Occlusal Guards