

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### 1.1. Product identifier

Cleanisept Wipes / Cleanisept Wipes Maxi

### 1.2. Relevant identified uses of the substance or mixture and uses advised against

#### Use of the substance/mixture

Disinfection and cleaning of medical inventory

### 1.3. Details of the supplier of the safety data sheet

Company name: Dr. Schumacher GmbH  
Street: Am Roggenfeld 3  
Place: 34323 Malsfeld / DEUTSCHLAND  
Telephone: +49 (0) 5664/9496-0  
Telefax: +49 (0) 5664/8444  
e-mail: post@schumacher-online.com  
Internet: www.schumacher-online.com

Responsible for the safety data sheet: sds@gbk-ingelheim.de

**1.4. Emergency telephone number:** INTERNATIONAL: +49 - (0) 6132 - 84463, GBK GmbH (24h - 7d/w - 365d/a)  
England and Wales: NHS Direct - 0845 4647; Scotland: NHS 24 - 08454 24 24  
24

## SECTION 2: Hazards identification

### 2.1. Classification of the substance or mixture

This mixture is not classified as hazardous according to Directive 1999/45/EC.

#### Classification according to Regulation (EC) No. 1272/2008 [CLP]

This mixture is not classified as hazardous according to Regulation (EC) No. 1272/2008.

### 2.2. Label elements

#### Additional advice on labelling

The product does not require a hazard warning label in accordance with EC directives/the relevant national laws.

### 2.3. Other hazards

Not known.

## SECTION 3: Composition/information on ingredients

### 3.2. Mixtures

#### Chemical characterization

Disinfectant tissues

## SECTION 4: First aid measures

### 4.1. Description of first aid measures

#### General information

Remove contaminated soaked clothing immediately.

If you feel unwell, seek medical advice.

#### After inhalation

Move to fresh air in case of accidental inhalation of vapours or decomposition products.

In the event of symptoms refer for medical treatment.

#### After contact with skin

Wash off immediately with soap and plenty of water.

Consult a doctor if skin irritation persists.

#### After contact with eyes

Rinse immediately with plenty of water, also under the eyelids.

If eye irritation persists, consult a specialist.

#### After ingestion

Do not provoke vomiting. Consult physician. Attention in case of vomiting - acute danger of suffocating, produced by foaming ingredients. Rinse mouth. Make drink some glasses of water. The decision whether to provoke vomiting is to be taken by a physician.

#### **4.2. Most important symptoms and effects, both acute and delayed**

Contact with eyes, skin or oral tissues may cause irritation.

#### **4.3. Indication of any immediate medical attention and special treatment needed**

Treat symptoms.

---

### SECTION 5: Firefighting measures

#### **5.1. Extinguishing media**

##### **Suitable extinguishing media**

Product does not burn, fire-extinguishing activities according to surrounding.

##### **Unsuitable extinguishing media**

Full water jet.

#### **5.2. Special hazards arising from the substance or mixture**

Fire may produce:

Carbon monoxide (CO), carbon dioxide (CO<sub>2</sub>) and nitrogen oxides (NO<sub>x</sub>).

Chlorine compounds.

#### **5.3. Advice for firefighters**

Use breathing apparatus with independent air supply.

Protective suit.

##### **Additional information**

Cool containers at risk with water spray jet.

Fire residues and contaminated firefighting water must be disposed of in accordance with the local regulations.

---

### SECTION 6: Accidental release measures

#### **6.1. Personal precautions, protective equipment and emergency procedures**

In case of vapour formation use respirator.

Ensure adequate ventilation.

Use personal protective clothing.

#### **6.2. Environmental precautions**

Do not discharge into the drains/surface waters/ground water.

#### **6.3. Methods and material for containment and cleaning up**

Take up mechanically and collect in suitable container for disposal.

#### **6.4. Reference to other sections**

Observe protective instructions (see Sections 7 and 8).

Information for disposal see section 13.

---

### SECTION 7: Handling and storage

#### **7.1. Precautions for safe handling**

##### **Advice on safe handling**

Keep container tightly closed.

Use only in thoroughly ventilated areas.

Avoid contact with skin, eyes and clothing.

##### **Advice on protection against fire and explosion**

No special protective measures against fire required.

#### **7.2. Conditions for safe storage, including any incompatibilities**

##### **Requirements for storage rooms and vessels**

Keep containers tightly closed in a dry, well-ventilated place.

Keep at temperatures between 5°C and 40°C.

##### **Advice on storage compatibility**

Incompatible with oxidizing agents.

#### Further information on storage conditions

Keep away from food, drink and animal feeding stuffs.

#### 7.3. Specific end use(s)

Disinfection and cleaning of medical inventory

### SECTION 8: Exposure controls/personal protection

#### 8.1. Control parameters

#### 8.2. Exposure controls

##### Appropriate engineering controls

Ensure adequate ventilation, especially in confined areas.

##### Protective and hygiene measures

Wash hands before breaks and immediately after handling the product.

When using do not eat, drink or smoke.

Avoid contact with skin, eyes and clothing.

Remove and wash contaminated clothes before re-use.

##### Eye/face protection

Protective goggles in the event of risk of splashes.

##### Hand protection

Suitable are gloves made of following materials (breakthrough time  $\geq$  8 hours):

Natural rubber/latex - NR (0,5 mm) (use powder free and allergene-free products)

Polychloropren - CR (0,5 mm)

Nitrile rubber/nitrile latex - NBR (0,35 mm)

Butyl rubber - Butyl (0,5 mm)

Fluoro-rubber - FKM (0,4 mm)

Polyvinyl chloride - PVC (0,5 mm)

##### Skin protection

Long sleeved clothing (EN 368).

##### Respiratory protection

No personal respiratory protective equipment normally required.

### SECTION 9: Physical and chemical properties

#### 9.1. Information on basic physical and chemical properties

Physical state:	Non-woven fabric, impregnated	
Colour:	White	
Odour:	Mild	
pH-Value:	6 - 8	*)
Melting point:	< - 10 °C	*)
Initial boiling point and boiling range:	approx. 100 °C	*)
Flash point:	n.a.	
Lower explosion limits:	n.a.	
Density (at 20 °C):	approx. 0,99 g/cm <sup>3</sup>	*)
Water solubility: (at 20 °C)	Miscible	
Ignition temperature:	n.a.	
Viscosity / dynamic:	approx. 9 mPa·s	*)

#### 9.2. Other information

\*) All data refer to the solution

### SECTION 10: Stability and reactivity

#### 10.1. Reactivity

No decomposition if stored and applied as directed.

## **10.2. Chemical stability**

Stable under normal conditions.

## **10.3. Possibility of hazardous reactions**

Reactions with oxidizing agents.

## **10.4. Conditions to avoid**

To avoid thermal decomposition, do not overheat.

## **10.5. Incompatible materials**

Strong oxidizing agents.

## **10.6. Hazardous decomposition products**

Carbon monoxide (CO), carbon dioxide (CO<sub>2</sub>) and nitrogen oxides (NO<sub>x</sub>).

Chlorine compounds.

---

# **SECTION 11: Toxicological information**

## **11.1. Information on toxicological effects**

### **Acute toxicity**

Based on available data, the classification criteria are not met.

No toxicological data available.

### **Irritation and corrosivity**

Based on available data, the classification criteria are not met.

### **Sensitising effects**

Based on available data, the classification criteria are not met.

### **STOT-single exposure**

Based on available data, the classification criteria are not met.

### **Severe effects after repeated or prolonged exposure**

Based on available data, the classification criteria are not met.

### **Carcinogenic/mutagenic/toxic effects for reproduction**

Based on available data, the classification criteria are not met.

Carcinogenicity: Not classified.

Mutagenicity: Not classified.

Reproductive toxicity: Not classified.

### **Aspiration hazard**

Based on available data, the classification criteria are not met.

### **Additional information on tests**

Classification in compliance with the assessment procedure specified in the Regulation (EC) no 1272/2008.

## **Practical experience**

### **Other observations**

Contact with eyes, skin or oral tissues may cause irritation.

If appropriately handled and if in accordance with the general hygienic rules, no damages to health have become known.

---

# **SECTION 12: Ecological information**

## **12.1. Toxicity**

Ecological data are not available.

## **12.2. Persistence and degradability**

The surfactant(s) contained in this preparation complies(comply) with the biodegradability criteria as laid down in Regulation (EC) No.648/2004 on detergents. Data to support this assertion are held at the disposal of the competent authorities of the Member States and will be made available to them, at their direct request or at the request of a detergent manufacturer.

## **12.3. Bioaccumulative potential**

No data available.

## **12.4. Mobility in soil**

No data available.

## **12.5. Results of PBT and vPvB assessment**

No data available.

## **12.6. Other adverse effects**

Low hazard to waters.

### **Further information**

Do not flush into surface water or sanitary sewer system.

## **SECTION 13: Disposal considerations**

### **13.1. Waste treatment methods**

#### **Advice on disposal**

Where possible recycling is preferred to disposal.

Can be incinerated, when in compliance with local regulations.

#### **Waste disposal number of waste from residues/unused products**

070699 WASTES FROM ORGANIC CHEMICAL PROCESSES; wastes from the MFSU of fats, grease, soaps, detergents, disinfectants and cosmetics; wastes not otherwise specified

#### **Contaminated packaging**

Empty containers should be taken for local recycling, recovery or waste disposal.

Contaminated packaging should be emptied as far as possible and after appropriate cleansing may be taken for reuse.

Packaging that cannot be cleaned should be disposed of like the product.

## **SECTION 14: Transport information**

**Land transport (ADR/RID); Marine transport (IMDG); Air transport (ICAO); Inland waterways transport (ADN)**

### **14.1. UN number:**

No hazardous material as defined by the transport regulations.

### **14.2. UN proper shipping name:**

No hazardous material as defined by the transport regulations.

### **14.3. Transport hazard class(es):**

No hazardous material as defined by the transport regulations.

### **14.4. Packing group:**

No hazardous material as defined by the transport regulations.

### **14.5. Environmental hazards**

No hazardous material as defined by the transport regulations.

### **14.6. Special precautions for user**

No hazardous material as defined by the transport regulations.

### **14.7. Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code**

No hazardous material as defined by the transport regulations.

## **SECTION 15: Regulatory information**

### **15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture**

#### **EU regulatory information**

1999/13/EC (VOC): 0 %

#### **National regulatory information**

### **15.2. Chemical safety assessment**

For this substance a chemical safety assessment has not been carried out.

## **SECTION 16: Other information**

### **Abbreviations and acronyms**

ADR = Accord européen relatif au transport international des marchandises Dangereuses par Route

RID = Règlement concernant le transport international ferroviaire de marchandises dangereuses

ADN = Accord européen relatif au transport international des marchandises dangereuses par voie de navigation intérieure

IMDG = International Maritime Code for Dangerous Goods

IATA/ICAO = International Air Transport Association / International Civil Aviation Organization

MARPOL = International Convention for the Prevention of Pollution from Ships

IBC-Code = International Code for the Construction and Equipment of Ships Carrying Dangerous Chemicals in Bulk

GHS = Globally Harmonized System of Classification and Labelling of Chemicals

REACH = Registration, Evaluation, Authorization and Restriction of Chemicals

CAS = Chemical Abstract Service

EN = European norm

ISO = International Organization for Standardization

DIN = Deutsche Industrie Norm

PBT = Persistent Bioaccumulative and Toxic

LD = Lethal dose

LC = Lethal concentration

EC = Effect concentration

IC = Median immobilisation concentration or median inhibitory concentration

#### **Further Information**

Data of items 4 to 8, as well as 10 to 12, do partly not refer to the use and the regular employing of the product (in this sense consult information on use and on product), but to liberation of major amounts in case of accidents and irregularities.

The information describes exclusively the safety requirements for the product(s) and is based on the present level of our knowledge.

The delivery specifications are contained in the corresponding product sheet.

This data does not constitute a guarantee for the characteristics of the product(s) as defined by the legal warranty regulations.

(n.a. = not applicable; n.d. = not determined)

---

*(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*