



Dryz® Blu post rinse
PHOTO: Dr. F. Shull, SC, USA.

Create a Dry Working Environment Every Time

Dryz[®] blu

NEW

Gingival Hemostatic Retraction Paste



- New 'blu' color makes it even easier to see the material.
- Stops gingival bleeding and seepage quickly.
- Aluminum Chloride acts fast and won't slough tissue like Ferric Sulfates.
- New creamy consistency adapts well and washes off clean.
- Available in syringes and single dose capsules, which eliminate waste and keep asepsis.

New Dryz Blu features a bright, contrasting aqua blue color making the material more visible against tooth structure. The slight stiffness of the paste provides better displacement of tissue when the material is pressed subgingivally (with either a compression cap or cotton roll).

The material rinses easily and because of the aqua blue color, any residual is quickly spotted and removed simply with a water/air spray.



3.9

SKU #	Product
S186	NEW Dryz® Blu Intro Kit Includes (30) unit dose capsules (0.25 g/capsule) & Dispensing Gun.
S187	NEW Dryz® Blu Unit Dose Capsules Includes (30) unit dose capsules (0.24 g/capsule).
S190	NEW Dryz® Blu Syringes Includes (7) 0.5ml syringes and (15) applicator tips.
S191	NEW Dryz® Blu Syringes Value Pack Includes (25) 0.5ml syringes and (55) applicator tips.
S181	Dryz® Luer Lock Applicator Tips Includes (40) applicator tips.